

I.V. House® UltraDressing® 330 Series

Protection for Neonates and Infants



I.V. House UltraDressing 330S



I.V. House UltraDressing 330M



I.V. House UltraDressing 330L

The I.V. House UltraDressing 330 Series is designed to protect the IV insertion sites of the tiniest patients. The clear, plastic, ventilated dome^{20,22} paired with the comfortable fabric wrap provides IV insertion site protection for very low birth weight neonates up to normal birth weight neonates. The 330 Series can be used on the hand, arm, foot, or leg.

Innovative Design

The I.V. House UltraDressing 330 Series covers the IV insertion site to prevent movement of the catheter.²⁰ This is important because movement can lead to complications, including accidental dislodgement, which can cause serious delay of treatment.

It is the position of the Infusion Nurses Society (INS)¹⁹ that when an infusion is running (continuous or intermittent) peripheral IV insertion sites should be routinely assessed for redness, tenderness, swelling, drainage, and or/paresthesia, numbness, or tingling at the specific frequency depending on age of patient. The use of the I.V. House UltraDressing 330 Series permits easy visualization of the IV insertion site allowing for a more frequent check without disturbing the patient.

I.V. House UltraDressing 330 Features

Adjustable, resealable Velcro® tab eliminates over-taping of fragile skin.

Smooth edges and a soft, pliable fabric wrap for safe, gentle IV protection.

Easy visibility and monitoring of IV insertion site to prevent injuries such as infiltration and extravasation.^{2, 3, 5, 6, 13, 14, 15, 18, 19, 20, 21, 22, 25, 26, 27, 33, 34, 35}

Eliminates exposed looping, a primary cause of dislodgement.²⁰

Maximizes dwell time, minimizes the need for painful, traumatic restarts.^{8, 16, 23, 31}

Offered in three sizes, the design fits either hand, arm, foot, or leg.



I.V. HOUSE®
Protection Over and Above™

Visit www.ivhouse.com for Instructional Videos, Directions, Product Details, and Ordering Information. Or call: **800-530-0400**

© 2015 I.V. House, Inc. All rights reserved. *Footnotes available at ivhouse.com/references

I.V. House® UltraDressing 330 Series > Product Specifications

