

Your saline flush is now your first line of defense against contamination.


SwabFlush®

Offers clinicians the convenience of having SwabCap ready when they need it, built into their saline flush!

Use SwabFlush for all your saline flush needs.

When used between medications, the SwabCap may be removed and saved for later use.

Simple to use!



As a final flush, SwabFlush ensures that the disinfection cap is always readily available to protect the connector between line accesses.



SwabCap
LUER ACCESS
VALVE CAP WITH 70% IPA
AS A DISINFECTANT

SwabCap® will disinfect the valve five (5) minutes after application and maintains a disinfected valve surface for up to seven (7) days if not removed!

Single Use Only
Sterile Packaging
Luer Lock Design

PRODUCT DESCRIPTION

SwabFlush® 10/10mL Saline Flush Syringe
SwabFlush® 5/10mL Saline Flush Syringe
SwabFlush® 3/10mL Saline Flush Syringe

Code

E0103-01
E0103-55
E0103-35

Units/Box

100
100
100

Units/Case

600
600
800



Excelsior Medical, LLC
a Medline Company

Call toll free 800-487-4276 or visit www.excelsiormedical.com to find out more about SwabCap®

SwabCap®


Valve disinfection is a critical part of IV catheter maintenance and infection prevention.

Let SwabCap® help you:

- Achieve disinfection and protection of needleless connectors
- Standardize disinfection technique and drive compliance to protocol
- Visually monitor and measure compliance
- Save nursing time

**LUER ACCESS VALVE CAP
WITH 70% IPA AS A
DISINFECTANT**

Single Use Only
Sterile Packaging
Luer Lock Design



SwabCap. . . Designed to meet your Highest Quality Standards

Sterile Packaging

No cross contamination fears for your patient

Thread Cover Seal

Alcohol stays in

Medical Grade Sponge Pad


Does not break down over time in alcohol

Soft and Low Profile

Gentle on patient's skin

Always at your fingertips

Integrated into your flush syringe



Patent pending **thread cover design**
Alcohol stays in and keeps contaminants out

SwabFlush®

The **best** of both technologies



**EXCELSIOR
MEDICAL**
Excelsior Medical, LLC
a Medline Company

Excelsior Medical, LLC
Corporate Headquarters
1933 Heck Avenue, Neptune, NJ 07753
www.excelsiormedical.com

Toll Free 800-487-4276
Main Phone 732-776-7525
Fax 732-776-7600
E-mail: info@excelsiormedical.com

SwabFlush® is a registered trademark of Excelsior Medical, LLC a Medline Company.
SwabCap® is a registered trademark of ICU Medical, Inc.

SFLUSH sell sheet U.S. Rev 07